

## Biological Indicator Incubator

### B80



#### Description:

The B80 Biological Indicator Incubator is compact, versatile, reliable and easy to use. It is designed for single temperature incubation of the ethylene oxide Biological Indicators. The thermometer is included with the incubator, providing easy verification of incubation temperature.

#### Features:

- ▶ Energy-efficient design.
- ▶ Applicable to 24V power adapter.
- ▶ Accurate to  $\pm 0.3^{\circ}\text{C}$ .
- ▶ Process of experiment is visible.
- ▶ Electronic calibration. User can reset to any local standard.
- ▶ Power supply provided is UL, CSA and CE approved.
- ▶ Custom blocks made to fit any sample holder.

#### Technical Specifications:

Model	B80
Temp. control range	RT+5°C~100°C
Timing range	1min~99h
Display accuracy	0.1°C
Heating time	$\leq 6$ min(from 20°C~56°C)
Dimension(mm)	145×113×110
Temp. accuracy	$\leq \pm 0.3^{\circ}\text{C}$
Temp. uniformity	$\leq \pm 0.3^{\circ}\text{C}$
Temp. stability	$\pm 0.3^{\circ}\text{C}$
Voltage	DC24V/1.67A
Weight(kg)	0.5



## Optional blocks:

Type	Diameter	Capacity
MD40	Φ10mm	15
MD41	Φ12mm	12
MD42	Φ15mm	8
MD43	Φ16mm	8
MD45	Φ9mm	18

