## DZG SCIENCE

Mini Vortex Mixer

**MIX1000** 

## **Description:**

MIX1000 Mini Vortex Mixer is designed for "touch" mixing a variety of test tubes. With a 4.5mm orbit and adjustable speed of 1000-2800rpm, the mini vortex mixer instantly vortexes even the largest samples, including nearly full 50ml tubes. The small footprint allows it to be used on even the most crowded bench.

## **Features:**

- ▶ Shaking speed: 1000-3500rpm.
- > Powerful and reliable motor designed for daily laboratory use.
- ▶ "Touch" pressure activated operation.
- ► Accept a variety of test tubes, including the tubes up to 50ml.
- $\blacktriangleright$  Mini-sized, allowing it to fit on even the most crowded bench

## **Technical Specifications:**

Model	MIX1000
Speed	1000-3500rpm
Type of motor activation	Touch pressure activated operation
Type of movement	Orbital
Orbit	4.5mm
Max test tube size	Ф30mm
Max load	0.1kg
Voltage	AC100-240V, 50/60Hz
Power	24W
Dimension(mm)	125×115×76
Weight(kg)	0.55

