

Vortex Mixer

MIX1001



Features:

- ▶ Small and compact appearance, easy to carry;
- ▶ Brushless DC motor, strong power and maintenance-free;
- ▶ Imported silicone vibrating head, strong toughness, sufficient flexibility and long service life;
- ▶ The sample can be swirled quickly within 5 seconds;
- ▶ Press to start, stepless speed regulation, convenient operation and wide application;
- ▶ Using flexible anti-vibration foot pads, metal base, and triangular structure, the instrument has good stability, low operating noise and no displacement during operation.

Technical Specifications:

Model	MIX1001
Speed Range	500-3500rpm
Drive Mode	Brushless DC motor
Operation Mode	Press and jog operation
Circumference diameter	4.5mm
Tube diameter	Max.35mm
Gross Weight	780g
Net Weight	540g
Maximum processing capacity	50ml
Power	5W
Voltage	DC12V 1A
Dimension	110x10x75mm
Package Dimension	210x155x135mm

