

## Shaker for Microplates

### MS100



#### Description:

MS100 Shaker for Microplates is a type of high speed shaker, which is with brushless DC motor and microprocessor controlled technique, specially for elisa plates(96/384 wells), tissue culture plates(24/48/96 wells).

#### Features:

- ▶ LCD display system status and parameters.
- ▶ Stable and reliable shaking.
- ▶ Easy to operate with one touch knob.
- ▶ Automatic fault detection and buzzer alarm function.

#### Technical Specifications:

Model	MS100
Timing range	1min ~ 99h59min
200~1500rpm	200~1500rpm
Mixing orbit	3mm
Speed step value	10rpm
Sample capacity	4plaka
Power	45W
Voltage	AC100V~240V, 50/60Hz
Dimension(mm)	280×270×110
Weight(kg)	4.0

